Deutsch

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English

## Instructions For Use

This information is for the correct use of **Mindable: Social Phobia**.

Please read this information carefully and follow the instructions contained therein.

If any questions or problems arise, the Mindable team can be reached at the contact details below.

### 1. Intended Purpose

**Mindable: Social Phobia** is a digital health application and is designed for individuals from 18 years old who suffer from the symptoms of social phobia.

The application aims to help these individuals cope with their disorder and alleviate their symptoms.

The application is designed to offer and teach these individuals appropriate disorder-related content, methods, and exercises. These are based on established psychotherapeutic guidelines, approaches, and procedures.

The application can be used independently by individuals as a supplement to their usual health care.

The application is not a substitute for treatment by a health care provider (such as a primary care physician, specialist, or psychotherapist) and is not intended to provide information for diagnostic or treatment decisions.

### 2. Minimum Age

The application is meant for affected individuals aged 18 years and older.

### 3. Application Area

The application is designed to support people diagnosed with symptoms of social phobia (ICD-10: F40.1).

### 4. Mode and Duration of Use

#### Program Content & Functionality

The program implemented in the digital application is based primarily on recognized treatment methods of cognitive behavioral therapy (CBT).

It is based on the treatment manual "Social Anxiety Disorder" (Stangier et al., 2016) and "Understanding and Changing Social Anxiety" (Hoyer, Härtling, & Sarnowsky, 2019) and was developed in close collaboration with the authors Prof. Dr. Stangier and Prof. Dr. Hoyer.

The application comprises several modules that build on each other.

First, users are provided with relevant information about the disorder and the cognitive explanatory models. With the help of questionnaires and individual questions, which would also be used in a normal psychotherapy, users are guided to develop an individual disorder model (anxiety model), to question dysfunctional expectations, as well as to identify and reflect on maintaining mechanisms of action.

In the exercise area of the application, attention training and behavioral experiments are instructed for implementation. On the one hand, users should learn to direct their attention in a targeted manner by means of these exercises, and on the other hand, users should continuously practice the reduction of safety behavior and self-attention, and question cognitive fears.

The goal is to provide users with a tool to help them manage their fears.

Furthermore, there is the possibility of self-assessment of the course of the sufferings with the help of targeted surveys. These are logged and summarized for the user in time-dependent graphs.

For support, if desired, users are encouraged by reminders and notifications to actively engage with the course content, to perform the exercises regularly and to adhere to set exercise goals.

Users can export recorded data in human- and machine-readable format.

Personal data is processed in accordance with the EU General Data Protection Regulation. Further information on data protection can be found in the data privacy policy.

The application can be used completely independently to bridge the waiting time for treatment, therapy-accompanying, as well as for relapse prevention.

### **Use Frequency & Duration**

To achieve long-term success, we recommend using the application intensively for at least 12 weeks. Thereby program content of the application can be worked at your own pace.

During the period of use, we recommend using the application at least once a week and in the best case daily.

Moreover, prolonged use is recommended so that behavioral changes and lifestyle adaptations can be solidified.

### **Conditions of Use**

The application can be used via a functioning smartphone or tablet with Internet access. An overview of the supported operating systems can be found at: [www.mindable.health/impressum](http://www.mindable.health/impressum)

After successful installation of the mobile application from the app store, an account must be created to use the application.

A valid e-mail address is required to create an account. This does not have to include a clear name and identity.

An internet connection is required during the app use. It is recommended to use only trustworthy networks.

The application is available in three languages (German, English and French). Good knowledge of one of the available languages is assumed.

### **Clinical Benefit**

As a medical device, the application can contribute to the reduction of the severity of social phobia symptoms.

## **5. Safety Information**

### **Contraindications**

Under what conditions should the application not be used?

The application is not allowed to be used in case of suicidal risk or acute self-injurious tendencies.

#### Explanation:

*Suicidal risk is present if there was a serious intention to take one's own life in the past 90 days (regardless of whether there was a specific plan to do so).*

### **Safety Advice & Side Effects**

It may happen that not every user can achieve symptom reduction with the application's program. This can potentially cause a sense of disappointment.

Dealing with upsetting issues can be stressful. If the burden of using the program is felt to be too great, the application should be paused immediately, the use of the program should be suspended for some time and the treating physician / psychotherapist should be contacted.

Not every exercise is equally suitable for every user. If an exercise causes difficulties or considered to be impairing, this exercise should not be continued, and the treating physician / psychotherapist should be contacted.

If any of the following occur during use, the treating physician / psychotherapist should be consulted immediately:

- if suicidal thoughts or thoughts of harming oneself arise or occur acutely

- if there is no improvement in the user's mental state over a longer period or if the mental state worsens acutely

If any unknown side effects occur related to the use of the application, please contact your physician / psychotherapist or the Mindable team immediately at: [support@mindable.health](mailto:support@mindable.health) (see further contact options below).

### Interactions

When narcotic substances such as benzodiazepines (e.g., alprazolam, clonazepam, diazepam, lorazepam), alcohol, or drugs are used simultaneously with the performance of behavioral experiments, the effects may cancel each other out. The central depressant effect of these substances opposes the conscious experience of the fear reaction in exercises and behavioral experiments. It is therefore not recommended to use the application under the influence of such substances.

#### Notice

Any change in prescribed medication and/or treatment should not be made without consulting the physician / psychotherapist in charge.

Please contact the physician / psychotherapist in charge of your treatment to find out about possible interactions.

### Emergency Information

In case of suicidal threat, acute crises, or emergencies, contact the physician / psychotherapist in charge immediately or dial the emergency numbers applicable in the country:

Germany:

- Emergency Call / Service: 112
- Police: 110
- Suicide / Mental Health Service: 0800 111 0 111, 0800 111 0 222

Austria:

- Emergency Call / Service: 144
- Police: 133

- Suicide / Mental Health Service: 142

Europe-wide emergency call: 112

For more information on country-specific emergency numbers and addresses, visit [adducation.info](http://adducation.info) and [suicide.org](http://suicide.org)

## 6. Further Information

### Contact Information

For technical support, please email:

[support@mindable.health](mailto:support@mindable.health)

Further information and contact options can be found on our website ([www.mindable.health](http://www.mindable.health)).

This instructions for use can be provided free of charge in printed form upon request to the above-mentioned support.

Should serious incidents occur in the context of application use, please contact Mindable Health GmbH immediately. In addition, the competent authority BfArM must be informed, which will be done by Mindable Health GmbH as the manufacturer.

### Software Status

These instructions for use refer to software version v2.x.x. The current software version of Mindable: Social Phobia can be found in the settings or in the Imprint: [www.mindable.health/impressum](http://www.mindable.health/impressum)

## Manufacturer

Mindable Health GmbH  
 Neue Grünstraße 17  
 D-10179 Berlin  
[www.mindable.health](http://www.mindable.health)

## Symbol Explanation

	Manufacturer
	Consult instruction for use
	Caution
	Medical Device
	Tradename
	Unique Device Identifier
	European conformity marking

## Status of Information

REV002-03-2023\_EN

## CE Marking of Conformity

Mindable: Sozial Phobia is a CE-marked medical device according to EU regulation 2017/745 (EU-MDR).

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informée, ce qui sera fait par Mindable Health GmbH, en tant que fabricant.